

# **EU**Vision

**COMPACT 22"**  
FULL HD Veterinary Endoscope System



## **New Features**

With **video recording playback function**

Professional **medical grade monitor and drivers**


## **Product Features**

- 1/2.3 SONY CMOS FULL HD Image Sensor
- Can be matched with stroboscopic light source for stroboscopic laryngoscopy
- FHD camera head and LED cold light source integrated
- Built-in USB FULL HD image storage, Integrated photo/video recording system
- Universal C-mount threading
- IPX8 waterproof camera

# EUVision

## COMPACT 22"

FULL HD Veterinary Endoscope System



### 22" Medical HD monitor

Extra anti-glare monitor Screen protector glass, 98% light transmittance

### LED medical light source

80W, support full range dimming, lift time  $\geq 50000$  Hours

### Camera head

Full HD imaging (1920 x 1080), 1/2.3" SONY CMOS, IPX8 waterproof, Full HD optical coupler C-mount standard

### Technical Parameters

Resolution: 1920(H) x 1080 (V) 60 frames

Screen size: 22"

Image format: JPG 1920 x 1080

Video format: Full HD 1920 x 1080

Light source: 80W LED cold light source

Product size: 528(L)x105(W)x370(H)

Weight: 9.85KG

**EurekaVet**<sup>®</sup>  
**MEDICAL**

Via 4 Novembre, 18 28040 PARUZZARO (NO)  
Tel. 0322 538109 info@eurekavet.it www.eurekavet.it